

A RANGE OF PRODUCTS THAT FACILITATE MORE EFFECTIVE SURGERY




A commitment to the highest quality and an ability to exceed international standards is what makes Frankenman and our Chex range of products stand out. We put health first. This manifests itself in our total dedication to a zero defect outcome. Our manufacturing, production, research and development, and customer service facilities are the envy of the international medical device sector.

We source and use only the highest quality materials and components from around the world. This means that not only do we adhere to internationally recognized measures of quality, such as CE and ISO 13485 certification, we regularly exceed them. Our commitment to the ongoing training and development

of our staff is second to none and our design capabilities are of the highest standards throughout the industry culminating in the award of full ISO 9001: 2000 certification for our Quality Management systems.



PROFESSOR ULRICH BAUMGARTNER OF THE KREISKLINIK OTTOBUREN PRESENTS HIS EXPERIENCE OF THE USE OF THE CHEX CPH RANGE OF STAPLERS



PROFESSOR ULRICH BAUMGARTNER

CHIEF OF THE DEPARTMENT OF GENERAL, VISCERAL SURGERY AND COLO-PROCTOLOGY SPECIALIST FOR SURGERY AND VISCERAL SURGERY, PROCTOLOGY (SPECIALIST AREAS OF INTEREST: COLORECTAL SURGERY)

Professor Baumgarnter is a recognised and experienced trainer in Stapled Hemorrhoidopexy, STARR and other colorectal procedures and has produced an excellent video guide on the techniques he employs.

WHEN YOU PUT PATIENTS' HEALTH FIRST, REMARKABLE THINGS FOLLOW



State-of-the-art manufacturing facilities.



OUR RANGE OF PRODUCTS INCLUDE

SINGLE USE RELOADABLE LINEAR CUTTER STAPLER



SINGLE USE RELOADABLE LINEAR STAPLER



SINGLE USE CURVED INTRALUMINAL CIRCULAR STAPLER



INFORMED DESIGN. SCRUPULOUSLY APPLIED. METICULOUSLY PRODUCED.

FOR MORE INFORMATION OR TO VIEW OUR ENTIRE PRODUCT RANGE PLEASE VISIT OUR WEBSITE.



SUITE B 13/F WING TAT COMMERCIAL BUILDING 121-125 WING LOK STREET SHEUNG WAN HONG KONG
T +852 3106 3035 F +852 3585 0519
E info@frankenman.hk

www.frankenman.hk

Chex and Chex Healthcare are registered trademarks of Frankenman International Limited



INFORMED DESIGN. SCRUPULOUSLY APPLIED. METICULOUSLY PRODUCED.



CIRCULAR STAPLERS FOR RECTAL PROLAPSE AND HEMORRHOIDS



A RANGE OF PRODUCTS TO MEET THE NEEDS OF MODERN SURGICAL PRACTICE

The introduction of the Chex™ CPH 34 and CPH 32 range of Circular Staplers for Stapled Hemorrhoidopexy has been received with the highest praise from Surgeons around the world. Undoubtedly...

“suma cum laude”

CPH 34 ACTUAL SIZE

21.60cm³ INTERNAL HOUSING VOLUME

TOP VIEW



1

SMOOTH, ERGONOMIC "LIP FREE" CONICAL DESIGN OF ANVIL.

2

LONGER "ANVIL ROD" OF ONE PIECE CONSTRUCTION ENABLES BETTER ACCESS AND VISUALIZATION WHEN SECURING "PURSE STRING" AND GREATER MECHANICAL STABILITY FACILITATING ACCURATE STAPLE FORMATION

3

FIRST TO INTRODUCE A MINIMUM TISSUE GAP OF 0.75mm AS STANDARD

4

CPH 34 FEATURES 4 SUTURE CHANNELS PROVIDING MORE OPTIONS FOR ACCURATE AND TENSION FREE PLACEMENT OF FULL SUTURES

10

"WING NUT" TO ASSIST CONTROL AND DEXTERITY WHEN OPENING AND CLOSING ANVIL

CPH 32 ACTUAL SIZE

17.17cm³ INTERNAL HOUSING VOLUME

SIDE VIEW



1

2

5

ONLY PRODUCT THAT HAS 32 STAPLES IN 2 ROWS TO ASSIST MORE EFFECTIVE HEMOSTASIS

6

HIGHEST QUALITY SURGICAL STEEL CUTTING BLADE


7

UNIQUE "CUTTING" MECHANISM WITH CIRCUMFERENTIAL PUSH ASSEMBLY

8

RADICAL RE-DESIGN OF STAPLE HOUSING GIVING SURGEON ACCURATE CONTROL OF TISSUE DRAWN INTO HOUSING

LOWERED TISSUE STOP TO FACILITATE SMOOTH CONTROLLABLE ENTRY OF TISSUE INTO SUTURE HOUSING



5

6

7

8

ACCESSORIES



9

FIRST TO SUPPLY TWO ANAL DILATORS. IN ADDITION TO STANDARD ANAL DILATOR ALL PRODUCTS ARE SUPPLIED WITH "BUTTERFLY" DILATOR ASSISTING WHEN ACCESS IS RESTRICTED

STAPLE INFORMATION

MODEL	CARTRIDGE COLOUR	STAPLE ROWS	OUTSIDE DIAMETER	CUTTING DIAMETER	STAPLE LEG LENGTH	STAPLE CROWN LENGTH	STAPLE COUNT	PER BOX
CPH 34	GREEN	2	34.5mm	26mm	4.2mm	3.8mm	32	3
CPH 32	WHITE	2	32.5mm	24mm	4.2mm	3.8mm	32	3

KEY FEATURES

- LARGE INTERNAL HOUSING CAPACITY TO ACCOMMODATE MORE TISSUE
- ADDITIONAL STAPLES TO PROVIDE EFFECTIVE HEMOSTASIS
- WING NUT CLOSING HANDLE FACILITATES EASE OF CLOSURE

A RANGE OF PRODUCTS THAT, QUITE SIMPLY, ARE NOW CONSIDERED THE BEST

When you put health first, you test everything as strenuously as possible. You test it yourself and have others test it for you, you ask questions, you go back if necessary, you incorporate what you've learned into the manufacturing process.

You devote all your energies to product quality and end-user cost reduction. Only then can you be satisfied that the products you produce are truly the best possible.

At Frankenman, we believe we're not so much in the business of making instruments but in the business of making surgeries more successful.

That's why we listen to doctors who've put our instruments to use. We read their reports. We invite them to think with us, to brainstorm, to help us develop new, innovative products.

That way we know we will be able to deliver the highest quality to the medical profession both today and way into the future.

