

When you put patients' health first, remarkable things follow.




**CHEX**<sup>™</sup>  
HEALTH FIRST

## SINGLE USE RELOADABLE LINEAR STAPLER

*"The staples are closer together, overlap more and are made from marginally thicker titanium wire. These factors suggest the CHEx stapler can produce a more intimate joint between two wound edges."*

University of Dundee, Comparative Assessment



## SINGLE USE RELOADABLE LINEAR STAPLER

INFORMED DESIGN.  
SCRUPULOUSLY APPLIED.  
METICULOUSLY PRODUCED.

- KEY FEATURES**
- + LARGE PARALLEL JAW OPENING: 20mm  
• Consistent compression • Consistent staple formation
  - + MANUAL/AUTOMATIC TISSUE LOCATING PIN
  - + RELOADS FOR DIFFERENT TISSUE THICKNESS
  - + CARTRIDGE SAFETY INTERLOCK DEVICE  
Prevents firing of used/empty cartridges
  - + MIDDLE INTERLOCK POSITION

### STAPLER SPECIFICATION

PRODUCT CODE				CLOSED STAPLE SIZE	STAPLING LINE LENGTH	STAPLE QUANTITY
LINEAR STAPLER	CARTRIDGE RELOADS	LINEAR STAPLER	CARTRIDGE RELOADS			
LS 30W	WHITE	LSR 30W	WHITE	1.0mm	32mm	23
LS 30G	GREEN	LSR 30G	GREEN	2.0mm	32mm	13
LS 30B	BLUE	LSR 30B	BLUE	1.5mm	32mm	13
LS 45G	GREEN	LSR 45G	GREEN	2.0mm	46mm	19
LS 45B	BLUE	LSR 45B	BLUE	1.5mm	46mm	19
LS 60G	GREEN	LSR 60G	GREEN	2.0mm	60mm	25
LS 60B	BLUE	LSR 60B	BLUE	1.5mm	60mm	25
LS 90G	GREEN	LSR 90G	GREEN	2.0mm	90mm	37
LS 90B	BLUE	LSR 90B	BLUE	1.5mm	90mm	37

**Frankenman International Limited,**  
Suite B, 13/F, Wing Tat Commercial Building,  
121-125 Wing Lok Street, Sheung Wan, Hong Kong  
t:+852 3106 3055 f:+852 3585 0519  
e: info@frankenman.hk

Chex and Chex Healthcare are registered trademarks of Frankenman International Limited.



www.frankenman.hk

**CHEX**<sup>™</sup>  
HEALTH FIRST

